

您的健康 我们呵护! Your Health We Care!



FLEXIBLE-Y ™ TYPE NEGATIVE SUCTION SHEATH

Flexible-Y [™] Type Negative Suction Sheath



SCOPE OF APPLICATION

This product is used in endoscopic urological surgery to help dilate the ureter to build access for insertion and removal of the endoscope and associated surgical instruments as well as fluid injection.

PRODUCT FEATURES

Passive bending

 Flexible-Y[™] front end facilitates to passively bend via endoscope to observe the stone in narrow renal calyx and improve visual field.



Flexible Sheath

Hydrophilic coating

 Hydrophilic coating at both inner and outer tube to avoid damage to the urinary tract and facilitate the discharge of calculus fragments.

High efficiency

 Clear the stone during smashing to save the surgery time, meanwhile, improve the stone clearance rate.

Security

 Rigid sheath, effectively reduce intra-operative flushing pressure and improve safety.

(Soft and smooth design)

 Flexible tip and smooth transition of connection port to protect ureter and device from damage during access.

Y Type					W Type				
Item	ID	OD	Length (after bending)	Length	ltem	ID	OD	Length (after bending)	Length
ZSR-UAY-0935	9 Fr	11 Fr	35CM	39CM	ZSR-UAW-1035	10 Fr	12 Fr	35CM	39CM
ZSR-UAY-0945	9 Fr	11 Fr	45CM	49CM	ZSR-UAW-1045	10 Fr	12 Fr	45CM	49CM
ZSR-UAY-1035	10 Fr	12 Fr	35CM	39CM	ZSR-UAW-1135	11 Fr	13 Fr	35CM	39CM
ZSR-UAY-1045	10 Fr	12 Fr	45CM	49CM	ZSR-UAW-1145	11 Fr	13 Fr	45CM	49CM
ZSR-UAY-1135	11 Fr	13 Fr	35CM	39CM	ZSR-UAW-1235	12 Fr	14 Fr	35CM	39CM
ZSR-UAY-1145	11 Fr	13 Fr	45CM	49CM	ZSR-UAW-1245	12 Fr	14 Fr	45CM	49CM
ZSR-UAY-1235	12 Fr	14 Fr	35CM	39CM	ZSR-UAW-1335	13 Fr	15 Fr	35CM	39CM
ZSR-UAY-1245	12 Fr	14 Fr	45CM	49CM	ZSR-UAW-1345	13 Fr	15 Fr	45CM	49CM
ZSR-UAY-1335	13 Fr	15 Fr	35CM	39CM	ZSR-UAW-1435	14 Fr	16 Fr	35CM	39CM
ZSR-UAY-1345	13 Fr	15 Fr	45CM	49CM	ZSR-UAW-1445	14 Fr	16 Fr	45CM	49CM
ZSR-UAY-1435	14 Fr	16 Fr	35CM	39CM					
ZSR-UAY-1445	14 Fr	16 Fr	45CM	49CM					



Dongguan ZSR Biomedical Technology Company Limited

Add: Unit 448, Qingfeng Road, Taihu Village, Sanzhong Village Committee, Qingxi Town, 523651 DONGGUAN City Guangdong Province P.R. China Fax: +86 769-87328976 Postcode: 523651

Email:sales@zsrbiomedical.com

www.zsrbiomedical.com

