


**DISPOSABLE
ENDOSCOPE STONE EXTRACTOR**

Disposable Endoscope Stone Extractor



SCOPE OF APPLICATION

During the operation, it is operated under endoscope and used to remove tissue, guide material and smash calculus

PRODUCT FEATURES

Multi-specification available

- Various models, three and four filaments design



Ergonomic handle

- Ergonomic handle design for easy operation

Isometric rotation


- Unique isometric rotation for easy stone removal

Nitinol basket


- Unique memorability and elasticity, anti-kink, easy to remove

Tipless design

- Tipless design to avoid trauma, reduce bleeding, and maintain a clear surgical field of vision.



Unique memory and flexibility



Tipless design

Multi-layer sheath

- Multi-layer sheath with smooth inner and outer wall to provide a smooth work channel, and with good balance and softness.

Model differentiation	Item	OD	Working Length
Three filaments	ZSR-SE01-1.7	1.7Fr	120cm
	ZSR-SE01-1.9	1.9Fr	120cm
Four filaments	ZSR-SE02-1.3	1.3Fr	120cm
	ZSR-SE02-1.9	1.9Fr	120cm
	ZSR-SE02-2.2	2.2Fr	120cm



Dongguan ZSR Biomedical Technology Company Limited

Add : Unit 448, Qingfeng Road, Taihu Village, Sanzhong Village Committee, Qingxi Town , 523651 DONGGUAN City Guangdong Province P.R. China
 Fax : +86 769-87328976 Postcode : 523651
 Email:sales@zsrbiomedical.com

www.zsrbiomedical.com

